

Have you seen this
revolutionary new
reading device
for
partially sighted
people?

VisioBook

VisioBook

Extremely small

Dimensions:

- folded: 340 × 310 × 54 mm
(appr. 13.4 (W) × 12.5 (D) × 2.1 (H) in.)
- in operation: 400 × 310 × 480 mm
(appr. 15.7 (W) × 12.5 (D) × 18.9 (H) in.)

VisioBook

Extremely light

Weight:

- 3,6 kg (appr. 7.9 lbs.) incl. power adapter

VisioBook

Extremely versatile

The VisioBook runs 5 hours on its internal battery!

You can use your VisioBook:

- in your home environment in different rooms
- in the classroom, as a student or a teacher
- in a library
- in the train or plane
- in your job on the desk or in meeting rooms

and you can stow it easily!

VisioBook

Extremely good!

a new high quality reading device suitable for almost any partially sighted user – from students to the active senior citizens!



- The first in a new class of portable magnifiers!
- Extremely lightweight, ultra-portable and easy to carry
- Easy setting up and folding down within a few seconds
- Aesthetically pleasing and space-saving design
- Intuitive operation and buttons with clear markings
- HD camera for both close range and distance viewing
- 5 hours of battery life!



Technical specifications

Easy operation:

- Button panel with 6 function keys and control dial for adjusting the size
- 7 color combinations each also inverse, full natural colors
- automatic focus can be switched off and on as desired
- Smart folding mechanism with 3 safe articulated joints

Monitor:

- Full-color 12.5" (32 cm)
- 16:9 format
- TFT wide screen with LED-backlight technology
- 1280 × 720 pixel (720p) HD camera for both close range and distance views

Magnification range:

- 1,8 to 30 times magnification

Settings:

- Magnification, dynamic contrast adjustment and picture optimising control

Size and weight:

- folded: 340 × 310 × 54 mm (appr. 13.4 (W) × 12.5 (D) × 2.1 (H) in.)
- in operation: 400 × 310 × 480 mm (appr. 15.7 (W) × 12.5 (D) × 18.9 (H) in.)
- appr. 3,6 kg (appr. 7.9 lbs.) incl. power adapter

Battery and voltage:

- Built-in Lithium-ion battery
- Operating time 5 hrs.
- External power adapter / battery charger 100-240V AC 50-60 Hz
- Power cons. in operation: < 15 W / in standby: < 1 W

Connections:

- USB
- HDMI (in and out)



Designed and made in Germany to high quality standards by:

BAUM Retec AG

In der Au 22
D-69257 Wiesenbach
Germany

Tel.: +49 6223 4909-0
Fax: +49 6223 4909-397

E-Mail: info@baum.de
Internet: www.baum.de